

# Low Temperature Chamber MN-3







# At a glance

# Your benefits

- ✓ Made in Germany parts and construction of high quality for high durability.
- Energy efficiency refrigeration plant and lighting provide optimal energy efficiency.
- ✓ Non corrosive all metal materials are galvanized and durable plastic-coated.
- High level of standard fitting:
  - $\circ$  Each lamp bank is dimmable separately in 0,5% steps as a standard.
  - Chamber parameters can be controlled and programmed with an industry-standard touch panel.
  - Each chamber comes with LAN connectivity as a standard for remote programming, parameter control and data logging. This also enables remote diagnosis und -service.
- ✓ **Sophisticated construction** We build the chambers at site out of 6 single pieces if necessary.

# Technical

- ✓ **Overall dimension** 900 x 900 x 2.090 mm (D x W x H).
- ✓ Inner dimension 1,59 m<sup>2</sup> working area on 3 tiers, 25 cm growing height per tier, given an equal lamp bank- and shelf distribution.
- ✓ Temperature from -10°C (without light) resp. 0°C (with light) until +40°C, given a maximum temperature variance of ±1,0°C.
- Air conditioning energy-efficient refrigeration system with hot-gas bypass-control and RPM-controlled ventilation fans. Dual evaporator-architecture for best possible temperature equability even during defrost periods.

#### ✓ **Flexible lighting:** the lamp banks are dimmable in 0,5%-steps and can be fitted with

- True Daylight white-LED
- True Daylight dual white LED
- True Daylight PLUS white-LED
- True Daylight dual PLUS white-LED
- Fluorescent lamps with various light intensities
- Optional spectrum enrichment with red- and far-red-LEDs or LEDs pickable from a broad nanometer range.
- Multichannel LED-panels with the light colors blue, white, red and far-red or many other nanometer ranges to choose from.
- Intuitive and comfortable operation industry-standard 12"-touchscreen at the chamber or remotely via standard network connection.





# General

poly klima®, a young and innovative company, designs and builds custom-made climatic walk-in rooms and growth chambers for environmental simulation for various research fields at universities and institutes.


Our team of experts benefits from many years of experience in conception and manufacturing of climatic- and lighting solutions for plant biological research.

The MN-series growth chambers from poly klima® excel with their technical equipment and their wide temperature range.

# poly klima<sup>®</sup> Low Temperature Chamber MN-3

The poly klima<sup>®</sup> Low temperature chamber MN-3 was developed for applications that require lower temperatures. Due to its flexibility it is also good for material-testing. The specially designed evaporator-architecture provides best possible temperature-equability even during defrost periods.

Modell MN-3 is loaded with three lamp banks and three shelves, that can be moved vertically within its tier. The shelves consist of white coated grating and can be pulled out on drawer rails.





### Design

All metal parts used are galvanized and coated with white, reflective durable plastic. Therefore, corrosion is not possible. The inner compartment is made out of stainless steel.

The chamber walls, the floor, the top cover and the door are made of steel sheet metal which is polyurethane foamed without any thermal bridges. Placed on the chamber floor there is a steel trough with a condensate drain. The chamber door is lockable.

The shelves are made out of white coated steel plates and can be slid out safely without danger of tilting. Shelves and lamp banks can be moved vertically within their tier. A tier is defined by the specially designed air conduit plates that can provide a vertical airflow. For constructional reasons the air conduit plates cannot be moved vertically.

The chamber is being delivered fully assembled.

The chamber is placed on (braked) casters and can be moved easily.

### Air-conditioning

Air-cooled refrigeration system, extremely energy-efficient with hotand gas control **RPM-regulated** ventilation fans, which significant minimizes wind stress for plants and energy consumption of the chamber during night periods. A doubleevaporator provides best possible temperature equability even during periods. defrost In case one evaporator freezes over the other evaporator takes over and the first one will be defrosted. To keep the



level of icing on the registers as low as possible, the working evaporators change every two hours.


This model is equipped with a variable air conduit. It can easily be switched from horizontal to vertical.


*Horizontal*: The air inside is travelling horizontal over the shelf and is led back up to the evaporator. There it is climatized and led back to the back-wall air-channel.

*Vertical*: With an easy to use sliding plate mechanism the back-wall air-holes can be closed and the openings for the air-conduit plates can be opened at the same time fort o realize vertical air-flow. This helps reducing condensate on cell culture lids like pedri dishes. For constructional reasons the air conduit plates cannot vertically be moved









For a free airflow through the compressor ventilation grill there should be at least 20 cm free space left/right to and

behind the chamber. We also are ready to adjust the ventilation grills placement to the local conditions. Furthermore, it is important that the chambers heat load and the heat load of other machines probably also placed in the room can be led away. As an alternative the room can be climatized, because ambient temperatures above 30°C will be critical to the chambers functionality and may result in a chamber shut-down.

Condensation water is being collected in the drain pan below the evaporator and led out of the chamber through the back wall. Also, a stainless-steel drain is placed on the chamber floor. From the integrated  $\frac{3}{4}$ " hose fitting there it can be led away to a floor drain or a condensate pump.

#### Temperature

Standard temperature range:  $-10/0^{\circ}C$  (without/with\* light) up to  $+40^{\circ}C$  given a maximum temperature variance of  $\pm 1,0^{\circ}C$ . (\* min. possible temperature withlight depending on build-in light intensity)

#### Dehumidification

Dehumidification with the evaporator in conjunction with an additional heating as a standard for rel. humidity values from ambient down to 45% r.H. ( $\pm$  5%; rel. humidity values valid in an ambient temperature range from 15°C bis 30°C; depending from ambient air humidity and light level inside the chamber).

### Lighting

The thorough and equidistant arrangement of the lamps on the light fixtures in conjunction with the white plastic coating inside the compartment ensures excellent light homogeneity over the whole growing area. There is no "fall-off" in the margin areas.

You have the choice between several lighting solutions: various white-LED-systems, fluorescentlamps or multichannel-LED solutions.

The spectrum of all below shown white-LED- and fluorescent-solutions can be enriched with further LEDs, selectable from a broad range of nanometers here (280 nm - 830 nm), like separate dimmable red-LED (e.g. 660nm) or far-red-LEDs (e.g. 730 nm).





#### White-LED-Solutions:

- True Daylight Standard White-LED
- White-LEDs with 3.000 K color-temperature.
- · Up to 400  $\mu$ mol/m<sup>2</sup>/s<sup>-1</sup> intensity (measured at 15 cm distance)
- Dimmable in 0,5%-steps from 100% to 1%.
- Color rendering index CRI=94.
- 120° radiation angle per LED.
- Very harmonic spectral response, ideal for plants like Arabidopsis.
- True Daylight dual White-LED
- White-LEDs with 2.700 K and 6.500 K color-temperature.
- · Each light color controllable separately.
- Up to 1.000  $\mu$ mol/m<sup>2</sup>/s<sup>-1</sup> intensity (measured at 15 cm distance)
- Dimmable in 0,5%-steps from 100% to 1%.
- · Color rendering index CRI=95.
- $\cdot$  120° radiation angle per LED.
- $\cdot$  Very harmonic spectral response, ideal for applications that require a higher light intensity or variable light quality.
- True Daylight PLUS White-LED
- White-LEDs with 4.000 K color-temperature.
- Up to 1.200  $\mu$ mol/m<sup>2</sup>/s<sup>-1</sup> intensity (measured at 15 cm distance)
- Dimmable in 0,5%-steps from 100% to 1%.
- Newest LED-technology
- · Color rendering index CRI=96,7.
- $\cdot$  120° radiation angle per LED.
- · Very harmonic and full spectral response, ideal for many applications.
- True Daylight dual PLUS white-LED
- Two types of white-LEDs with 3.000 K and 6.500K color-temperature.
- Each light color controllable separately.
- Up to 1.200  $\mu$ mol/m<sup>2</sup>/s<sup>-1</sup> intensity (measured in 15 cm distance)
- Dimmable in 0,5%-steps from 100% to 1%.
- Newest LED-technology
- · Color rendering index CRI=96,6.
- $\cdot$  120° radiation angle per LED.
- · Very harmonic and full spectral response, ideal for many applications, variable light-quality.

#### Multichannel-LED-Solutions

For plant-biological applications that require the control of different spectral parts or for other special applications, we offer our polyphoLED multichannel-LED Panels. Up to 12 different color channels can be addressed and controlled separately. There is a broad palette of narrow-banded LEDs available from 285nm up to 830nm. The resulting spectral response and light intensity are depending on the print configuration.








#### Fluorescent tubes

TL-D fluorescent tubes, light color 840 (neutral white) with good performance and energy efficiency over the chamber's temperature range.

The lamp banks are dimmable in 0,5% steps as a standard.



There are more

intensity groups to choose from (intensities measured at 15 cm distance):

a.) 5 to 250  $\mu$ mol/m<sup>2</sup>/s b.) 8 to 400  $\mu$ mol/m<sup>2</sup>/s

Of course, light intensities are customizable according customer wishes!

### Operation

Chamber parameters can be controlled and programmed with a high-grade industry-standard touch panel on the chamber or a remote computer. The 12" graphic display ensures a quick and intuitive programming and shows all actual and nominal values.

With the visualizing software on every accessauthorized Windows-based computer in the network all parameters can be displayed, edited and programmed comfortably from your desk, just as you would be standing in front of the chamber.

Temperature Conf Temperature SetPoint Temperatur	27.3 °C	Lighting   warm white 100.0 %   (cold white 50.0 %	
		(warm white 50.0 %) (cold white 50.0 %)	
Humidity Control Humidity SetPoint Humidity	88.5 % 84.7 %	(warm white 50.0 %) (cold white 50.0 %)	
Program Controller	CO2	Fresh air	Disinfection
rogram Controller ate : inactive ogram name :			

All alarm messages will be shown in text-messages

on the touch panel and additionally be forwarded via email or SMS.





### Options

- Ultrasonic humidification for humidity levels inside the chamber up to 85% r.H.\* (±5%, depending on ambience humidity and light level inside the chamber). \*The figures given for the relative humidity control are valid in an ambient temperature range from +15°C to +30°C. The humidification system will be automatically disabled at chamber temperatures ≤0°C.
- > <u>Reservoir humidification</u> for Entomology, to avoid interfering of ultrasonic on reared insects.
- > <u>Gas application</u> for chamber inside with  $CO_2$  or  $O_2$ .
- > <u>Entomology package:</u> chamber specially designed for rearing insects.
- > <u>Door-window</u> with or without light-tight cover.

This is just an extract of the most important options. Our chambers are practically customizable to every requirement. We look forward your challenge!

## Contact

We appreciate your interest and your questions! Just give us a phone call or drop us an email.

We are glad to offer advice and help at any time!

#### poly klima GmbH

Parkstraße 14 D-86462 Langweid-Foret • Germany Tel. +49 (0)821/ 650 752 75 • Fax +49 (0)8161/68 66 066

info@polyklima.de • service@polyklima.de



### www.polyklima.com

© poly klima 2019

